Part number: A20855  
Serial number:

# Universal Light Source

## User's Manual and Technical Reference

Part I: English  
Part II: German

122 Charlton Street, Southbridge, MA 01550-1960, USA  
Phone: 508-765-9744 Fax: 508-765-1299  
[www.us.schott.com/lightingimaging](http://www.us.schott.com/lightingimaging)  
Email: [lightingimaging@us.schott.com](mailto:lightingimaging@us.schott.com)

# SCHOTT

## **Warning**

1. Never look directly at the lamp when on.
  2. Lamp and surrounding surfaces may be hot!
  3. Do not use this unit near water or in an area with excessive moisture.
  4. Do not place flammable materials near the unit.
  5. Refer to the Introduction Section for safety instructions.
- 

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## 1. Introduction:

### 1.1. Product Description:

Congratulations on the purchase of your Universal Light Source. This light source provides cool, intense, white light for all types of fiber optic illumination components. The Universal Light Source is CSA listed to UL-1571 and CE compliant to low voltage and EMC directives. The Universal Light Source was designed and manufactured under the ISO 9000 guidelines.

The Universal Light Source was designed to be an extension of the DCR<sup>®</sup> III Light Source. However, any power supply that doesn't exceed 21 VDC at 8 amps will operate the Universal Light Source. The removable Modulamp<sup>®</sup> assembly contains the lamp holder, which supports three different lamp choices depending on intensity, lifetime and uniformity requirements, or any MR-16 lamp that does not exceeding the power requirements.

### 1.2. Safety Instructions:

1. Read, understand and follow all instructions in this manual.
2. Keep all safety and operating instructions for future reference.
3. Do not block ventilation openings on this unit. Do not impede airflow.
4. Use only standard glass type cleaners. Do not use solvents or petroleum distillates.
5. Never spill liquid on the unit or push objects into ventilation openings.
6. Disconnect power cord when unit will be left unused for long periods of time.
7. Allow the unit to cool before servicing.
8. Do not service the unit beyond what is described in this manual. Return the light source to an authorized SCHOTT service center. For a Return Material Authorization see Section 4.5.

## 2. Operation:


### 2.1. Voltage Input:

1. The Universal Light Source was designed to accept a varying DC voltage between 0 and 21 Volts. This voltage is lamp-dependent; therefore you should not exceed the lamp manufacturers rated voltage and power consumption.

**Note:** The Universal Light Source utilizes a DC fan that will vary with the voltage supplied to the unit.

## 2.2. Turning the Light Source “on”:

1. Insure the power switch on the back panel is in the “off” position (A).
2. Plug an appropriate cable into the Amp CPC Connector (B).
3. To help minimize current loss the Amp CPC Connector pins have been split into two groups: pins 1 and 2 are the +VDC (C), and Pins 3 and 4 are the – VDC (D).
4. All information regarding the AMP CPC Connector can be found on their web site: <http://www.amp.com>. Search under part number: 206430. Cabling options are available. Please contact SCHOTT Sales Department.

5. Make sure that the captive thumbscrew on the Modulamp<sup>®</sup> assembly is firmly tightened.
6. Insert any fiber optic component with a compatible light source adapter into the Modulamp<sup>®</sup> receptacle.
7. Tighten the fiber optic positioning thumbscrew (see diagram).
8. Turn the light source power switch to the “on” position.
9. To turn the unit off, press the power switch to the “off” position.

## 3. Modulamp<sup>®</sup>:

### 3.1. Modulamp<sup>®</sup> Assembly Removal:

1. Turn the power switch to the “off” position.
2. Remove the fiber optic component.
3. Allow the lamp to cool.
4. Unscrew the Modulamp<sup>®</sup> assembly captive thumbscrew until the head springs forward.

**CAUTION: Lamp, lamp socket and surrounding surfaces may be hot!**


5. Pull the top of the Modulamp<sup>®</sup> assembly forward and lift the unit out of the recess.
6. Reach into the recess and depress the release tab on the connector.
7. Slide the connector apart and remove the Modulamp<sup>®</sup> assembly.

### 3.2. Lamp Replacement:

1. Turn the power switch to the “off” position.
2. Remove the fiber optic component.
3. Allow the lamp to cool.
4. Unscrew the Modulamp<sup>®</sup> assembly captive thumbscrew until the head springs forward (H).



**CAUTION: Lamp, lamp socket and surrounding surfaces may be hot!**

5. Pull the top of the Modulamp<sup>®</sup> assembly forward and lift the unit out of the recess.
6. Push back on the lampholder release lever to raise the lamp (L).
7. Gently unplug the ceramic lamp socket (K).
8. Lift the release lever. Insert a new lamp into the lamp holder. Push the lamp all the way down. Be careful not to touch the bulb or the inside of the reflector. Fingerprints may affect the light output.
9. Push the ceramic socket on to the two lamp pins until flush with the lamp base.
10. Return the Modulamp<sup>®</sup> assembly into the recess and tighten captive thumbscrew.
11. For complete Modulamp<sup>®</sup> assembly removal and installation, see Sections 3.1 and 3.3.

### 3.3. Modulamp® Assembly Installation:

1. Make sure the power switch is in the “off” position.
2. Slide the connector into the mating receptacle until the release tab clicks into place.
3. Place the Modulamp® assembly into the recess with the two-bottom “ears” riding over the bottom edge of the recess first.
4. Slide the Modulamp® assembly into place.
5. Depress the captive thumbscrew and tighten firmly.

### 3.4. Modulamp® Options:

The Modulamp® assembly allows you the convenience of changing the lamp, lamp type or intensity control (iris) by simply exchanging Modulamp® assemblies.

<b>Modulamp Assembly</b>	<b>Part Number</b>
Without Iris	<b>A08301</b> , EKE Position <b>A08301.1</b> , EJA Position <b>A08301.2</b> , DDL Position
With Iris	<b>A08321</b> , EKE Position <b>A08321.1</b> , EJA Position <b>A08321.2</b> , DDL Position

### 3.5. Lamp Type Change:

1. Each Modulamp® assembly can be adjusted to accommodate any of the three recommended lamps (EJA, EKE or DDL). Repositioning after changing the lamp type is very important. It will assure maximum intensity at the fiber optic input.
2. After removing the Modulamp® assembly from the housing, note the two sets of pre-drilled holes provided in the Modulamp® base.
3. The set of holes closest to the front of the Modulamp® assembly provide the correct positioning for the EJA and EKE lamps. The set of holes furthest from the front of the Modulamp® assembly are positioned for the DDL lamp.
4. To move the lamp holder, remove the two screws holding the lamp holder to the Modulamp® base, align the set of holes for the desired position, and screw the holder back into place.

## 4. Appendix:

### 4.1. Cleaning:

The metal housing and the “Lexan” front and back decals have a durable finish that should retain their original luster for many years to come. Cleaning the exposed areas with a commercial glass cleaner will help to maintain the finish.

Press the power switch on the front of the unit to the “off” position. Remove the power cord from the CPC connector on the back of the unit. Wipe the exposed areas of the housing and front and back plates with a soft cloth or paper towel moistened with any commercial glass cleaner.

**CAUTION:** Do not use detergents, excess water, treated cloths, harsh cleaning agents or sprays. Use cleaning fluid sparingly. If fluid is spilled into the interior let the unit dry thoroughly before using.

Periodically, dust should be removed from the unit using commercially available cans of compressed air. Pay special attention to the fan blades at the rear of the unit, the air vents in the front and bottom of the unit and the Modulamp<sup>®</sup> assembly.

### 4.2. Troubleshooting:

If you are unsuccessful at resolving the following conditions, contact SCHOTT. See Section 4.5 for a **Return Materials Authorization (RMA)**. Do not attempt to repair the light source. Removing the cover will void the warranty.

Do not attempt to troubleshoot the light source with any remote control device connected. If the light source operates after removal of the remote connection, check your remote set-up.

1. Neither fan nor lamp operates:  
(Note that fan and lamp circuits are tied together.)
  - Make sure the power cord is inserted completely into the AMP connector and also into the correct power source.
  - Make sure AMP CPC is plugged in completely.
  - Check the extension cord for damage.
  - Check power switch.
2. Universal Light Source turns on/off:
  - Universal Light Source ambient too high.
  - Check vent holes for obstruction.

### 4.3. Product Specifications:

1. Physical Dimensions:  
Overall Height: 4.62" (177 mm)  
Overall Width: 7.27" (185 mm)  
Overall Depth: 8.61" (219 mm)  
Fiber Receptacle: .720" (18.3 mm) diameter
2. Environmental Conditions:  
Operating Temperature Range: 32 to 113° F, (0 to 45° C)  
Storage Temperature Range: -13 to 185° F, (-25 to 85° C)  
Relative Humidity Range: 0 to 95%, Non-condensing
3. Electrical Specifications:  
Input Voltage: 0 – 21 VDC @ 8 Amps

**CSA** listed to **UL** (1571)

**DECLARATION OF CONFORMITY**

APPLICATION OF COUNCIL DIRECTIVE(S): **73/23/EEC Low Voltage Directive**  
**89/336/EEC EMC Directive**

STANDARDS (S) TO WHICH CONFORMITY IS DECLARED: **EN61326-1: 1997/A1: 1998 Class A**  
**EN61000-3-2: 1995+A1: 1998+A2: 1998**  
**Product Safety: DCR®III, DCR®III Plus, DCR® IV,**  
**IEC 950, ACE® IEC 598-1**

TYPE OF EQUIPMENT: **Light Source, Halogen**

MODEL NO.: **ACE I, ACE REMOTE, DCR® III, DCR®III Plus DCR®IV,**  
**DCR®III Remote Light Source Kit, MHR®-50, Universal Light Source**

MANUFACTURER'S NAME: **SCHOTT North America, Inc.**

MANUFACTURER'S ADDRESS: **122 Charlton Street, Southbridge, MA 01550-1960**

SIGNATURE:  
NAME: **Heinz Kaiser** *Heinz Kaiser*  
POSITION: **V.P & General Manager**

IMPORTER'S NAME: \_\_\_\_\_  
IMPORTER'S ADDRESS: \_\_\_\_\_

**NOTE:** Improvements may result in specification or feature changes without notice.

#### 4.4. Replacement Parts List:

The following replacement parts are available from SCHOTT to maintain and support your Universal Light Source:

A08110, Lamp – EJA  
A08120, Lamp – EKE  
A08130, Lamp – DDL  
A29506, Lamp Socket , Assembly  
D08375.063, Scr, Assy, Thumb, Modulamp®  
D20855.14, Manual, Remote, Technical Reference.

**Please Note:** SCHOTT Light Sources are designed to be used exclusively with Ushio or GE lamps. Light Source performance will not be guaranteed with any other manufacturer's lamps.

#### **4.5. Support/Service/RMA Policy:**

Support: SCHOTT maintains a variety of support services to assist you:

**In USA: Telephone (1)315-255-2791, Fax (1)315-255-2695**  
**In Europe: Telephone +49-6131-667796, Fax +49-6131-667850**

Contact your local SCHOTT distributor from whom you made your purchase for an immediate response.

Be sure to have your part number and serial number available, as well as a complete description of the problem or situation for the quickest, most accurate assistance.

#### **Service/RMA Policy:**

Any service required for any reason must be performed by SCHOTT or an authorized service representative. All service outside warranty will be performed upon the purchaser's request according to normal service charges in effect at the time.

To return any item, whether for warranty repair or chargeable servicing, an RMA number (Return Materials Authorization) must be obtained from SCHOTT. This number must be clearly visible on the shipping label. All shipping must be prepaid.

SCHOTT guarantees all warranty repairs will be completed within two weeks of receipt. All units will ship prepaid using our shipping method of choice. Alternate shipping methods will be shipped freight collect.

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