



i-Scope® Viewer for Android™ fit for SCHOTT MORITEX skin counseling system

Software now available on Google Play™

Asaka (Japan), May 13, 2013 – SCHOTT MORITEX, a provider of cosmetic counseling systems, has registered its skin counseling system software i-Scope® Viewer for Android™ on Google Play™. The software can now be downloaded from Google Play™ for free, enabling users of i-Scope® USB 2.0 to easily observe and record skin data on Android™ tablet computers. SCHOTT MORITEX will exhibit i-Scope® USB 2.0 and i-Scope® Viewer for Android at CITE Japan 2013 (booth A7-2).

**Visit SCHOTT MORITEX:
Cosmetic Ingredients &
Technology Exhibition
Yokohama, Japan
May 15-17, 2013
Booth A7-2**

SCHOTT MORITEX skin counseling systems are used by international cosmetic manufacturers as marketing tools. The new i-Scope® Viewer for Android™ was developed to achieve a more handy counseling method using tablet computers. The product is now available on Google Play™ in three languages: English, Japanese and Chinese.

i-Scope® Viewer for Android™ is the fitting software for i-Scope® USB 2.0. The program allows users to easily observe their skin conditions in detail: by moving the scope along their skin, shooting images and saving them, comparing two images in one screen or viewing multiple saved images simultaneously. i-Scope® USB 2.0 is a scope equipped with a 1.3 mega CMOS sensor that can be used to observe skin, scalp and hair conditions by connecting the USB cable to a computer. Five lens options are available.

Whereas the software for i-Scope® USB 2.0 could only be used with Windows® before, it is now also compatible with Android OS



and can thus also be used on Android™ tablets. Users who had purchased i-Scope® USB 2.0 for Windows® in the past can now also download this software to their Android™ tablets.

Technical data:

Lenses:

- 30×PL lens: For skin observation of texture and pigmentation
- 50×PL lens: For skin observation of texture and pigmentation
- 120× lens: For scalp and hair observation
- 200× lens: For scalp and hair observation
- 700× lens: For hair and cuticle observation

Features of i-Scope® USB 2.0 Android™ and specification of i-Scope® Viewer for Android™

- Connection: Wired (mini-USB cables sold separately)
- No additional setting necessary on tablet PC; easy connection
- Stable view of moving pictures
- Compatible devices: Versions after Android ver. 3.2 are USB host capable (except for some devices)
- Can check skin conditions by shooting still images
- Can display two images in one screen. Can check different parts of the skin
- Can save displayed images
- Can display multiple saved images
- Available in English, Japanese and Chinese

** Windows is a registered trademark of Microsoft Corporation in the United States and other countries.*

** Android is a registered trademark of Google, Inc, registered in the U.S. and other countries.*

** NFC PLUG is a registered trademark of Dai Nippon Printing Co., Ltd.*

** i-Scope is a registered trademark of SCHOTT MORITEX Corporation. .*



Photo ID 202001: The software i-Scope® Viewer for Android™ can now be downloaded from Google Play™ for free, enabling users of i-Scope® USB 2.0 to easily observe and record skin data on Android™ tablet computers. Photo: SCHOTT MORITEX

Download link to these photographs in print quality:
<http://www.schott-pictures.net/presskit/201334.iscope>

More press photographs are available for downloading here:
www.schott-pictures.net

SCHOTT is an international technology group with more than 125 years of experience in the areas of specialty glasses and materials and advanced technologies. SCHOTT ranks number one in the world with many of its products. Its core markets are the household appliance, pharmaceuticals, electronics, optics, solar power, transportation and architecture industries. The company is strongly committed to contributing to its customers' success and making SCHOTT an important part of people's lives with its high-quality products and intelligent solutions. SCHOTT is committed to managing its business in a sustainable manner and supporting its employees, society and the environment. The SCHOTT Group maintains close proximity to its customers with manufacturing and sales units in 35 different countries. Its workforce of around 16,000 employees generated worldwide sales of approximately 2 billion euros for the 2011/2012 fiscal year. SCHOTT AG, with its headquarters in Mainz, Germany, is owned by the Carl Zeiss Foundation.

SCHOTT MORITEX, integrated into the SCHOTT Business Unit Lighting and Imaging, develops, manufactures and distributes lighting and imaging solutions and supplies optical components and systems for machine vision, digital imaging and medical applications. SCHOTT MORITEX Corporation is listed on the first section of the Tokyo Stock Exchange. SCHOTT AG holds 71.6% of the shares of SCHOTT MORITEX Corporation.



Media Contact:

SCHOTT AG
Corporate Public Relations
Dr. Haike Frank
Hattenbergstr. 10
55122 Mainz, Germany
Phone: +49 6131 66 4088
E-mail: haike.frank@schott.com
Internet: www.schott.com

Sales Contact:

Japan, Korea and Taiwan:

SCHOTT MORITEX Corporation
Cosmetic Instrument Sales, Medical, Regional Sales Japan / Korea
Masatomo Oikawa
3-13-45 Senzui, Asaka-shi, Saitama, 351-0024, Japan
Phone: +81-48-218-2538
Fax: +81-48-462-6714
E-mail: masatomo.oikawa@schott-moritex.com
Internet: www.schott-moritex.com

China:

Lighting and Imaging SCHOTT Shanghai
Amanda Zhang
Unit 301, Innov Tower, No.1801 Hong Mei Road
200233, Shanghai, China
Phone: +86- 21-3367-8000-146
Fax: +86- 21-3367-8080
E-mail: amanda.zhang@schott.com

Singapore:

Lighting and Imaging SCHOTT Singapore Pte. Ltd.
Hock-Boon Yee
8 Admiralty Street, #05-01 Admirax,
757438, Singapore, China
Phone: +65-6488-2366
Fax: +65-6286-0838
E-mail: sales.singapore@schott.com

Europe:

Lighting and Imaging SCHOTT France SAS
Maren Sowa
6 bis, rue Fournier
92110, FR-Clichy
Phone: +33-1 4087-3946
Fax: +33-1 4087-3990
E-mail: maren.sowa@schott.com



PRESS INFORMATION



USA:

SCHOTT North America Inc, San Jose, CA
Jason Baechler
6862 Santa Teresa Blvd.,
95119, San Jose
Phone: +1-408-363-2100
Fax: +1-408-363-9980
E-mail: cosmetics@us.schott.com